

SUBMED GOWN



SUBMED offers you wearable safety & protection solutions with SUBMED Gown series. SUBMED Gowns are designed based on international standards for all clinicians. Our production flexibility may also provide special requirements without compromising ergonomics, comfort and quality.

SUBMED Protection Material

Safe, Easy and Protection

SUBMED Surgical Gowns

SUBMED's gown fabrics offer high-level protection. Each gown is latex-free. We offer different fabric options for surgical gowns, each with a special style, function and reinforcement and meets or exceeds the minimum performance requirements of the EN 13795-1, EN 13795-2 standards.

The Gowns materials have a good hydrostatic head and the extra repellent treatment gives patients and healthcare professionals the protection they need.

High performance with high comfort.....!

All our surgical gowns offer operating theatre staff excellent comfort and ease of use. High-tech, breathable fabrics keep the wearer comfortable and dry, even during long procedures. The gowns feature sleeves with ultrasonic bonding, which contribute to greater strength and comfort.

Note :

BS EN 13795-2:2019
Surgical clothing and drapes. Requirements and test methods Clean air suits.

DIN EN 13795-1
Surgical clothing and drapes - Requirements and test methods - Part 1: Surgical drapes and gowns

Submed GOWNS

- SMMMS (35-40-43 gsm) Material with natural, fluid-repellent protection
- High performance with comfortable and breathable
- Sterile and non sterile versions available



STERILE REF NO	NON-STERILE REF NO	PRODUCT DESCRIPTIONS	STERILE UNITS PER CARTON	NON-STERILE UNITS PER CARTON
SG-001-1	NG-001-1	Standart Surgical Gown, S, 117,5 cm $\pm$ 5%	80 pcs	80 pcs
SG-001-2	NG-001-2	Standart Surgical Gown, M, 125 cm $\pm$ 5%	80 pcs	80 pcs
SG-001-3	NG-001-3	Standart Surgical Gown, L, 132,5 cm $\pm$ 5%	75 pcs	80 pcs
SG-001-4	NG-001-4	Standart Surgical Gown, XL, 140 cm $\pm$ 5%	70 pcs	80 pcs
SG-001-5	NG-001-5	Standart Surgical Gown, XXL, 150 cm $\pm$ 5%	65 pcs	80 pcs
SG-001-6	NG-001-6	Standart Surgical Gown, XXXL 155 cm $\pm$ 5%	60 pcs	80 pcs



CTN. DIMS.
600x400x500 mm

Submed GOWNS

- SMS/PE (35-40-43/28 gsm) Material with natural, fluid-repellent protection
- AAMI Level 3 protection
- High performance with comfortable and breathable
- Sterile and non sterile versions available



STERILE REF NO	NON-STERILE REF NO	PRODUCT NAME,SIZE	STERILE UNITS PER CARTON	NON-STERILE UNITS PER CARTON
SG-002-1	NG-002-1	Reinforced Surgical Gown, S, 117,5 cm $\pm$ 5%	80 pcs	80 pcs
SG-002-2	NG-002-2	Reinforced Surgical Gown, M, 125 cm $\pm$ 5%	80 pcs	80 pcs
SG-002-3	NG-002-3	Reinforced Surgical Gown, L, 132,5 cm $\pm$ 5%	75 pcs	80 pcs
SG-002-4	NG-002-4	Reinforced Surgical Gown, XL, 140 cm $\pm$ 5%	70 pcs	80 pcs
SG-002-5	NG-002-5	Reinforced Surgical Gown, XXL, 150 cm $\pm$ 5%	65 pcs	80 pcs
SG-002-6	NG-002-6	Reinforced Surgical Gown, XXXL, 155 cm $\pm$ 5%	60 pcs	80 pcs



CTN. DIMS.
600x400x500 mm

Submed GOWNS

- SMS/PE (35-40-43/28 gsm) Material with natural, fluid-repellent protection
- AAMI Level 4 protection
- High performance with comfortable and breathable
- Sterile and non sterile versions available



STERILE REF NO	NON-STERILE REF NO	PRODUCT NAME,SIZE	STERILE UNITS PER CARTON	NON-STERILE UNITS PER CARTON
SG-003-1	NG-003-1	Full Reinforced Surgical Gown, S, 117,5 cm $\pm$ 5%	80 pcs	80 pcs
SG-003-2	NG-003-2	Full Reinforced Surgical Gown, M, 125 cm $\pm$ 5%	80 pcs	80 pcs
SG-003-3	NG-003-3	Full Reinforced Surgical Gown, L, 132,5 cm $\pm$ 5%	75 pcs	80 pcs
SG-003-4	NG-003-4	Full Reinforced Surgical Gown, XL, 140 cm $\pm$ 5%	70 pcs	80 pcs
SG-003-5	NG-003-5	Full Reinforced Surgical Gown, XXL, 150 cm $\pm$ 5%	65 pcs	80 pcs
SG-003-6	NG-003-6	Full Reinforced Surgical Gown, XXXL, 155 cm $\pm$ 5%	60 pcs	80 pcs



Submed GOWNS

- SSMMS/PE (33/22 gsm) Material with natural, fluid-repellent protection
- AAMI Level 4 protection
- High performance with comfortable and breathable
- Sterile and non sterile versions available



STERILE REF NO	NON-STERILE REF NO	PRODUCT NAME,SIZE	STERILE UNITS PER CARTON	NON-STERILE UNITS PER CARTON
SG-004-1	NG-004-1	Extra Reinforced Surgical Gown, S	80 pcs	80 pcs
SG-004-2	NG-004-2	Extra Reinforced Surgical Gown, M	80 pcs	80 pcs
SG-004-3	NG-004-3	Extra Reinforced Surgical Gown, L	75 pcs	80 pcs
SG-004-4	NG-004-4	Extra Reinforced Surgical Gown, XL	70 pcs	80 pcs
SG-004-5	NG-004-5	Extra Reinforced Surgical Gown, XXL	65 pcs	80 pcs
SG-004-6	NG-004-6	Extra Reinforced Surgical Gown, XXXL	60 pcs	80 pcs



Submed GOWNS

- Spunbond (30-35-40 gsm) Material with natural, fluid-repellent protection
- High performance with comfortable and breathable
- Available in standart size



STERILE REF NO	NON-STERILE REF NO	PRODUCT NAME,SIZE	STERILE UNITS PER CARTON	NON-STERILE UNITS PER CARTON
-	NG-005-1	Patient Gown without Cuff, 30 gr/m ² , 120 cm ±5%	-	80 pcs
-	NG-005-2	Patient Gown without Cuff, 35 gr/m ² , 120 cm ±5%	-	80 pcs
-	NG-005-3	Patient Gown without Cuff, 40 gr/m ² , 120 cm ±5%	-	80 pcs
-	NG-006-2	Patient Gown with Cuff, 35 gr/m ² , 120 cm ±5%	-	80 pcs
-	NG-006-3	Patient Gown with Cuff, 40 gr/m ² , 120 cm ±5%	-	80 pcs



CTN. DIMS.
600x400x500 mm

Submed GOWNS

- Spunbond (22 gsm) Material with natural, fluid-repellent protection
- High performance with comfortable and breathable



STERILE REF NO	NON-STERILE REF NO	PRODUCT NAME,SIZE	STERILE UNITS PER CARTON	NON-STERILE UNITS PER CARTON
-	NG-009	Visitors Gown, 22 gr/m ² , 92 cm ±5%	-	180 pcs
-	-		-	
-	-		-	
-	-		-	
-	-		-	



CTN. DIMS.
600x400x500 mm

Submed GOWNS

- SMS (35-40-46 gsm) Material with natural, fluid-repellent protection
- High performance with comfortable and breathable
- 3 Pockets, V neck
- 2 Parts



STERILE REF NO	NON-STERILE REF NO	PRODUCT NAME,SIZE	STERILE UNITS PER CARTON	NON-STERILE UNITS PER CARTON
-	NSS-001-1	Scrub Suits, S	-	100 pcs
-	NSS-001-2	Scrub Suits, M	-	100 pcs
-	NSS-001-3	Scrub Suits, L	-	100 pcs
-	NSS-001-4	Scrub Suits, XL	-	100 pcs
-	NSS-001-5	Scrub Suits, XXL	-	100 pcs



Submed GOWNS

- Warm Up Jackets are used to reduce risk of infection and create a layer of protection against splash for surgeons in all kinds of operations
- Long and Short sleeves available
- High performance with comfortable and breathable
- 2 Pockets





NON-STERILE REF NO	PRODUCT NAME,SIZE	NON-STERILE UNITS PER CARTON
NWJ-004-1	Warm Up Jacket, S Long Sleeve Full Bottomed	110 pcs
NWJ-004-2	Warm Up Jacket, M Long Sleeve Full Bottomed	100 pcs
NWJ-004-3	Warm Up Jacket, L Long Sleeve Full Bottomed	90 pcs
NWJ-004-4	Warm Up Jacket, XL Long Sleeve Full Bottomed	80 pcs
NWJ-004-5	Warm Up Jacket, XXL Long Sleeve Full Bottomed	70 pcs



Common Features of Submed Surgical Gowns

According to area of use, material type Submed gowns can be:


- Ultrasonic bonded
- Stitched
- Double Stitched
- Adjustable Velcro
- Minimum 10 cm female
- Minimum 5 cm male
- Ultrasonic Bonded Belts




	STERILE & NONSTERILE		
	a	b	c
S	65	118 cm	150
M	66	125 cm	150
L	68	130 cm	150
XL	69	135 cm	150
XXL	70	140 cm	155
XXXL	71	145 cm	160

Product & Carton Labels


STERILE PRODUCT LABEL



NON-STERILE PRODUCT LABEL



STERILE CARTON LABEL



NON-STERILE CARTON LABEL

