

SUBMED MASK



SUBMED Protection Material

Safe, Easy and Protection

SUBMED Surgical Mask

SUBMED Masks are designed according to international Bacterial Filtration Efficiency levels. They ensure high comfort in breathing while reduce the transfer of particles. SUBMED Masks are available in various colour for adults and children.



Submed Masks

TYPE IIR Mask Specifications

Outer Material	Filter Layer	Inner Material	Nose-Piece	Earloop
Spunbond 25 gr	Melt Blown 25 gr	Spunbond 15 gr	Aluminium	Polyester and Polyurethane

Test	Type IIR
Bacterial filtration efficiency (BFE), % :	≥ 98
Differential pressure (Pa/cm²) :	< 60
Splash resistance pressure (kPa) :	≥ 16,0
Microbial cleanliness (cfu/g) :	≤ 30

Length : 175 mm
Width : 95 mm
Height of Earloop : 75 mm
Lenght of Earloop : 170 mm
Pleat Depth : 10 mm
Lenght of Nose-Piece : 85 mm
Color : Blue
Applied Standards : EN 14683+AC 2019 - TYPE IIR
Quality Certificates : EN ISO 13485:2016, EN ISO 9001:2015
Product Classification : Medical device acc. to MDD 93/42/EEC, CE class 1



Submed Masks

Surgical Face Mask TYPE IIR



Excellent
Breathability



Adaptable
Nose Bar



Perfect
Fitting



High Filtration
Capacity



Hypoallergenic

Packaging Types

Polyback of 10 pcs



Box of 50 pcs



REF	UNI./PCS./BOX	Ctn. Dims.
MS001-4	50 pcs/ box - 48 box/ ctn	60x40x27 cm
MS001-4(10)	10 pcs/ bag - 250 bag/ ctn	60x40x27 cm

TYPE II Kids Mask Specifications

Outer Material	Filter Layer	Inner Material	Nose-Piece	Earloop
Spunbond 25 gr	Melt Blown 25 gr	Spunbond 15 gr	Aluminium	Polyester and Polyurethane

Test	Type II
Bacterial filtration efficiency (BFE), % :	≥ 98
Differential pressure (Pa/cm²) :	< 40
Splash resistance pressure (kPa) :	≥ 16,0
Microbial cleanliness (cfu/g) :	≤ 30

Length :	145 mm
Width :	90 mm
Height of Earloop :	45 mm
Lenght of Earloop :	135 mm
Pleat Depth :	10 mm
Lenght of Nose-Piece :	85 mm
Color :	Blue
Applied Standards :	EN 14683+AC 2019 - TYPE II
Quality Certificates :	EN ISO 13485:2016, EN ISO 9001:2015
Product Classification :	Medical device acc. to MDD 93/42/EEC, CE class 1



Excellent
Breathability



Adaptable
Nose Bar



Perfect
Fitting



High Filtration
Capacity



Hypoallergenic

Surgical Face Msk TYPE II - Light Blue Color



Polyback of 10 pcs

Box of 50 pcs



REF	UNI./PCS./BOX	Ctn. Dims.
MS001-6	50 pcs/ box - 60 box/ctn	60x40x27 cm
MS001-6(10)	10 pcs/ bag - 250 bag/ctn	60x40x27 cm